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General instructions for Series SE 552 combination pH electrodes

Caution!
Please read these instructions carefully before putting the electrode into operation, in order to ensure faultless function. Operation and servicing of the electrode must only be performed by specially trained personnel authorized by the operating company.

1. Introduction
Knick pH electrodes of the SE 552 Series are low-maintenance electrodes with pressurized liquid electrolyte, silver-ion trap, and integrated temperature probe. They provide a temperature signal for automatic temperature compensation of the pH signal during calibration and operation. The SE 552 sensors are autoclavable and in-situ sterilizable.

2. Safety remarks
The electrodes are intended only for the use described in point 1 above. The electrodes consist partly of glass and should be handled carefully. If an electrode is to be cleaned or calibrated using acid or alkaline solutions, safety goggles and gloves should be worn. The electrodes may be operated and serviced only by personnel familiar with the equipment in question and who have read and understood these instructions. Please follow the Ex classification enclosed for an installation in a hazardous location. The electrodes are designed to withstand a maximum overpressure of 4 bars. Please take extra care during installation/fitting of the electrode to ensure that neither the polymer body/shaft of the electrode nor the thread (Pg 13.5) suffer any injury or damage.

3. Product description
The markings on each electrode or package label provide the following information:
Knick Electrode manufacturer
SE 552/x Electrode designation / Length coding
pH 0...14 pH measurement range
0...100 °C (140 °C) Temperature range for operation/sterilization
☞ II 1/2G EEx ia IIC T6/T5/T4/T3 Hazardous area classification
☞ SEV 05 ATEX 0143 X Designation and number of 0032 certificate
In addition, each electrode is provided with a serial number on the red part of the connector cap to allow identification / traceability.

4. Installation and preparation for use
1. On unpacking, check the electrode for possible mechanical damage of the pH sensitive glass, the junction, the shaft and the VarioPin (VP) connector. Report any damage immediately to your Knick service team.
2. Remove the watering cap and briefly rinse the electrode with de-ionized water. Carefully remove the protection droplet of silicone on the junction with the enclosed knife, as shown on the green pictures on the packing.
3. Again rinse the electrode well with de-ionized water. After rinsing, the electrode should only be dabbed dry with a tissue. Do not rub the pH-sensitive glass since this can lead to electrostatic charging and sluggish response times.
4. Carefully check the space behind the pH-sensitive glass for the presence of any air bubbles and remove them by gently swinging the electrode in a vertical plane.
5. Please refer to the instruction manual of the relative fitting in order to install the electrode correctly.
6. Connect the electrode to the pH transmitter, using a cable with VarioPin (VP) connector. Please refer to the wiring and terminal diagram accompanying the cable.

5. Operation
5.1 Calibration of the electrode and pH transmitter
For SE 552 electrodes, a 2-point calibration is recommended. After removing the watering cap, the electrode is dipped in succession into two different buffer solutions with given pH reference values and the pH transmitter is calibrated to these buffer values.

5.2 Sterilization of electrode
In the case of sterile processes, such as in fermentation applications, the electrodes are sterilized prior to the production cycle. Autoclaving of the whole fermenter/reactor unit including pH electrode is possible with this sensor type. Sterilization can also be effected in-situ by means of steam or superheated process medium in contact with the lower electrode parts.

5.3 Temperature probe
The integrated Pt 1000 temperature probe is primarily intended for automatic temperature compensation of the pH signal during calibration and normal operation, and not for the accurate indication or control of the process temperature.

6. Maintenance
1. Carefully rinse the electrode tip and junction with de-ionized water after each production cycle. Under no circumstances must measuring solution be allowed to dry on these parts!
2. When the electrode is not in operation, store it with electrode tip and junction well submerged in 3M KCl electrolyte.
3. If an electrode is stored dry for a few days by mistake, let it soak in electrolyte (3M KCl) for several hours before use.
4. Occasionally check the connector for possible traces of moisture. If necessary, clean well with de-ionized water or alcohol, and afterwards dry carefully.

7. Trouble-shooting
(sluggish response, increased resistance etc.)
At sluggish response time (e.g. after a long and dry storage period), reactivation of the pH sensitive glass may be necessary. Dip the pH sensitive glass (but no further electrode parts!) into reconditioning (HF) solution for 5 min. Afterwards, rinse well with de-ionized water and let electrode tip and junction soak in reference electrolyte for at least 12 hours. Rapidly increasing resistance of reference electrode and/or drift could be signs of protein precipitation. Use Pepsin/HCl cleaning solution to remove such precipitation.

8. Disposal
If no relative official regulations are in force, used or defective SE 552 electrodes and their packaging can be disposed of as ordinary household waste.

9. Warranty
In the event of manufacturing faults, 12 months warranty from date of delivery is granted.

English:
1) VP vonnector
2) Pg 13.5 thread
3) Washer, PTFE
4) O-ring (silicone), 10.5 x 2.7 mm
5) Electrode body
6) Junction
7) Glass membrane
8) Argenthal reference element
9) Liquid reference electrolyte
10) Silver-ion trap
11) Reference element
12) Temperature probe

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Elektronische Messgeräte GmbH & Co. KG